



USER GUIDE

DireCtQuant NAE

Nucleic acids electrophoresis reagent system

**2 X faster
run time**

(compared to TAE)

**2 X better
resolution**

(compared to TAE)

**2 X longer
usage time**

(compared to TAE)

**½ cheaper
per run**

(compared to TAE,
when used as indicated)

For Research Use Only

Not Intended for any Animal or Human Diagnostic or Therapeutic Use

Contents

ABOUT DireCtQuant NAE	2
COMPONENTS AND EQUIPMENT TO BE PROVIDED BY THE USER	3
GEL AND SAMPLE PREPARATION	4
GENERAL STATEMENT	6

ABOUT DireCtQuant NAE

DireCtQuant nucleic acid electrophoresis buffer (Ref #DCQNAE) is an innovative formula providing far superior speed and resolution comparing to established formulas as TAE, TBE, etc. With **DireCtQuant NAE** can be achieved twice as better resolution of DNA/RNA fragments and the separation time can be reduced at least two times using a very similar workflow and standard equipment. The buffer generates much less heat during the electrophoresis and is having very high buffering capacity which protects the sample during the electrophoresis. The buffer is supplied with 2X nucleic acid electrophoresis sample loading dye NAEL (Ref #DCQNAEL). This dye is having very low viscosity allowing easy sample preparation and loading and contains fast migrating red tracking dye .

One set of **DireCtQuant NAE** provides the necessary reagents for at least 250 complete mini gels (10 cm x 8 cm):

Ref. #DCQNAE 50X Electrophoresis buffer 1 x 2 L

Ref. #DCQNAEL 2X Sample loading buffer 1 X 20 ml

COMPONENTS AND EQUIPMENT TO BE PROVIDED BY THE USER

- Electrophoresis system capable of accommodating your gels

! IMPORTANT: Carefully check the maximum safety limits (voltage, current and power) of your electrophoretic equipment (power supply, electrophoresis chamber). Do not use this system with equipment that cannot safely provide the conditions listed in this protocol. Not following these precautions can result in poor performance and could lead to equipment failure and/or personal injury!

- Agarose
- Distilled Water
- Microwave or hot plate
- Pipette set
- Nucleic acid staining reagent
- Gel visualisation system
- Personal protective equipment

GEL AND SAMPLE PREPARATION

1. Dilute the electrophoresis buffer by diluting DireCtQuant NAE 50 times with distilled water. For example, for preparing 1 litre of 1X working buffer to 20 ml of 50X DireCtQuant NAE add 980 ml of distilled water. Mix well.
2. Use the prepared 1X buffer to prepare your agarose gel and to fill the electrophoresis tank. You can include nucleic staining reagent in your gel or stain after electrophoresis.
3. Add to your sample equal volume of DireCtQuant NAEL loading buffer 2X (e.g. for 20 µl sample add 20 µl 2X DCQNAEL). Mix well.
4. Load the prepared samples on the gel.
5. Apply up to 15 V/cm (for example using mini gel equipment with the distance between the electrodes of 15 cm, you can apply 225V).

! IMPORTANT: *Never exceed the maximal voltage/current limits of your equipment. This can result in damage of your equipment or personal injury!*

6. Follow the separation by the migration of the red front. The migration of the dye is indicative for the migration of 50 bp dsDNA. Stop the electrophoresis when the tracking dye reaches the end of the gel. Full separation is usually achieved in 10-15 minutes.
7. Visualize the nucleic acids bands as with gel documentation system of your choice. The loading buffer DCQNAEL contains unique fast, non-fluorescent migrating dye. It migrates in front of the smallest DNA/RNA fragment resolvable on agarose electrophoresis. It doesn't create any bands on the gel images, making possible correct detection and densitometry analysis of the whole gel. In contrast Bromphenol blue and Xylene cyanol FF based loading buffers exhibit 3-4Kbp. and 200-300bp. bands interfering with fragment detection and quantification.

- 8. Due to very high buffering capacity the NAE buffer can be reused at least 3X times without detectable changes in performance, pH and conductivity. Mix the buffer in the electrophoresis cuvette if planned to be reused.**

The system is compatible with all common nucleic acids dyes (ethidium bromide, SYBR Safe, GelRed*, *etc.*) and downstream applications (purification from gel, Northern Blot, Southern Blot, *etc.*)

**SYBR safe and GelRed are a trademarks of Invitrogen and Biotium respectively.*

GENERAL STATEMENT

Handle our products in accordance with safe laboratory practices:
wear suitable protective gloves, safety goggles and protective clothing.

Please read and understand the MSDS data provided and updated at directquant.eu

For Research Use Only

Not Intended for any Animal or Human Diagnostic or Therapeutic Use

In case of any doubt please feel free to contact DireCtQuant team at:
customersupport@directquant.eu